ENDOTINE® Midface-ST 4.5


REF CFD-060-0197

Instructions for Use


MICROAIRE®

For Surgery. For Life.™

1

2

3


7


8

4


5

6

9

10

Name	Ref # (ISO 7000) ²	Symbol	Definition	Use Standard
REF (Catalog #)	2493		<ul style="list-style-type: none"> Indicates the manufacturer's catalog number so that the medical device can be identified. Per EN980:2008, the REF symbol may be used without surrounding box. 	ISO 15223-1:2012 ¹
Registered Trademark	N/A		N/A	International standardized symbol
Sterilized using Irradiation (gamma)	2502		Indicates a medical device that has been sterilized using irradiation (gamma). Use of the Sterile symbol requires a use-by date (see Use-By Date symbol).	ISO 15223-1:2012 ¹
Quantity	N/A		Quantity	N/A
Lot / Batch Code	2492		Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1:2012 ¹
Caution	0434A /0434B		Indicates the need for the user to consult the Instructions For Use (IFU) for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the device itself.	ISO 15223-1:2012 ¹
Manufacturer	3082		<ul style="list-style-type: none"> Indicates the medical device manufacturer. This symbol shall be accompanied by the name and address of the manufacturer. The date of manufacture may be combined with this symbol. When using MicroAire as the manufacturer, use the MicroAire LLC symbol. 	ISO 15223-1:2012 ¹

Date of Manufacture	2497		<ul style="list-style-type: none"> Indicates the date when the medical device was manufactured. The date is expressed as YYYY-MM (e.g. 2015-11) or YYYY-MM-DD (e.g. 2015-11-29). If the symbol is filled (see Manufacturer symbol), both the date of manufacture and the name/address of the manufacturer may be combined in one symbol. 	ISO 15223-1:2012 ¹
Consult Instructions For Use (IFU)	1641		Indicates the need for the user to consult the Instructions For Use (IFU). Not required in conjunction with the Caution symbol, if applicable.	ISO 15223-1:2012 ¹
Do Not Reuse	1051		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure.	ISO 15223-1:2012 ¹
Temperature Limitation	0632		Indicates the temperature limits to which the medical device can be safely exposed. The upper and lower limits to temperature shall be indicated adjacent to the upper and lower horizontal lines.	ISO 15223-1:2012 ¹
Prescription	N/A	R_{only}	Caution: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician (or properly licensed practitioner).	FDA, Title 21, Chapter 1, Subchapter H, Part 801.15F
Do Not Use if Package is Damaged	2606		Indicates a medical device that should not be used if the package has been damaged or opened. This symbol may also mean "Do not use if the product sterile barrier system or its packaging is compromised".	ISO 15223-1:2012 ¹
Use-By Date	2607		Indicates the date after which the medical device is not to be used. This symbol shall be accompanied by a date to indicate that the medical device should not be used after the end of the month shown. The date is expressed as YYYY-MM (e.g. 2015-11) or YYYY-MM-DD (e.g. 2015-11-29).	ISO 15223-1:2012 ¹

¹ ISO 15223-1:2012 – "Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements"

² ISO 7000 – "Graphical symbols for use on equipment – Registered symbols"

ENDOTINE Midface-ST 4.5

Available Models

Part No.	Tine Length	Part Name
 CFD-060-0197	4.5 mm	ENDOTINE Midface-ST 4.5

DESCRIPTION

The ENDOTINE Midface-ST 4.5 consists of insertion tools and a bioabsorbable implant. The implant consists of a fixation platform attached to an anchoring leash (Figure 1). This device and its insertion tools are supplied sterile for single use only.

INDICATIONS/INTENDED USE

The ENDOTINE Midface™-ST 4.5 is intended for use in subperiosteal midface suspension surgery to fixate the cheek subdermis in an elevated position.

CONTRAINDICATIONS

1. Situations where internal fixation is otherwise contraindicated (e.g., infection).
2. Patients appearing to have very thin soft tissues of the face in which the implant may be visible or palpable.
3. Any known allergy or foreign-body sensitivities to plastic biomaterial.

PACKAGING, LABELING AND STERILIZATION

The ENDOTINE Midface-ST 4.5 and the supplied insertion tools are pre-sterilized by gamma irradiation and intended for single use only. Do NOT resterilize this device or the insertion tools. The device and its insertion tools should be accepted only if the factory packaging and labeling arrive intact. Do not use if packaging shows evidence of puncture(s), tampering, water contamination or other damage.

STORAGE INSTRUCTIONS

1. Store at room temperature $\leq 79^{\circ}\text{F}$ / 26°C in a dry place. Do not use beyond the expiration date listed on the label. Do not use if the temperature indicator dot has turned black; the indicator changes color if the device is subjected to temperatures that could damage it, even if the exposure is of short duration, as might occur during transportation.

INSTRUCTIONS FOR USE

1. The implant device (Figure 1) is labeled for orientation including the fixation platform (1A) and anchoring leash (1B).
2. The implant assembly (Figure 2) is labeled for orientation including the release mechanism (2A), anchoring leash (2B) and protective cover (2C).

3. The midface dissection from a temporal approach—with or without additional incisions (buccal, periorbital)—is performed using the surgeon's typical approach. In order to allow the fixation platform to lie over the maxillary antrum after elevation (instead of over the zygoma where it is likely to be visible), the dissection must continue to the inferior maxilla over the buccal sulcus.
4. Insert the implant assembly through the temporal incision into the desired position in the midface (Figure 3) with the protective cover oriented towards the skin as illustrated. If significant resistance to the assembly is encountered, do not force it into position. Remove the assembly and extend or enlarge your plane of dissection, with care to protect surrounding nerve structures, in order to provide adequate clearance for proper positioning of the platform at the inferior recess of the dissected tissue.
5. When comfortable with the fixation platform position at the inferior recess of the midface dissection, retract the protective cover by squeezing the release mechanism (Figure 4). Apply digital pressure over the cheek to engage the tines in the midface tissue (Figure 5), and remove the insertion tools (Figure 6).
6. Alternatively, the implant device may be removed from the insertion tools and inserted in a retrograde fashion through an oral incision. (Figure 7)
7. Elevate the midfacial tissues to the desired position by applying tension to the anchoring leash (Figure 8) exiting the temporal incision. At this point, if the fixation platform appears centered over the zygomatic prominence, undesirable skin irregularities or excessive palpability/sensitivity may result. Beginning the elevation with the device initially more inferior on the midface may avoid this problem.
8. Once the desired position is reached, exercise caution with excessive tissue manipulation or endoscopic visualization, as this may dislodge the fixation platform from the tissue.
9. After elevating the device and midface to the desired position, place suture through the leash (Figure 9) and into the temporal fascia where the soft tissue is of adequate thickness.
10. Trim the anchoring leash so that the end lies comfortably within the temporal incision. (Figure 10)
11. If adjustment or removal of the implant is desired, one may remove the suture and extract the fixation platform under direct endoscopic guidance. Alternatively, direct visualization of the platform is possible through the buccal incision, if present.
12. The patient should be instructed to exercise caution postoperatively for several weeks to avoid excessive manipulation or contact with the face.
13. Should removal of the device be required in the postoperative period, a small temporal incision to release the leash and an oral approach to remove the platform may offer the most straightforward method.

WARNING	Used to indicate that the safety of patients and hospital personnel could be involved.
CAUTION	Used to indicate procedures that must be followed to avoid damaging an instrument.
NOTE	Used to indicate the easiest means of carrying out the techniques.

WARNINGS

1. DO NOT RESTERILIZE the ENDOTINE Midface™-ST 4.5 or its insertion tools.
2. Do not use if storage temperature exceeds 110°F (43°C) or the temperature indicator on the outer box has turned black.
3. Discard open, unused devices.
4. This device and tools are for SINGLE USE ONLY and MUST NEVER BE RE-USED. Resterilization may cause loss of function or failure of the device. Risks of re-use may include patient infection due to lack of sterility.
5. Inspect all components for damage before use. Do not attempt to repair damaged components.

CAUTIONS

1. Inform the patient of possible adverse effects including palpability of the device prior to resorption, discomfort, visibility of the device under the skin, and surgical risks inherent to the procedure (to include nerve injury, facial musculature dysfunction, infection, scarring, device erosion through the skin, etc.). Some devices may require removal prior to their absorption due to discomfort, infection, reaction or other concerns.
2. Transient local fluid accumulation and/or sinus tract formation may occur. Aspiration may yield traces of implant material and may result in healing of the sinuses without adverse effects.

ADVERSE EFFECTS

Material sensitivity/allergic reactions in patients following surgery should be reported. Implantation of foreign materials in tissues can result in histological reactions.

Caution: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

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**MicroAire Surgical Instruments, LLC**