

www.merillife.com

Corporate Headquarters, Vapi

Meril Healthcare Pvt. Ltd.
H1-H3, Meril Park
Survey No. 135/2/B & 174/2
Muktanand Marg, Chala
Vapi 396 191. Gujarat. India.
T: +91 260 3052 100
F: +91 260 3052 125
E: orthopedics@merillife.com

International Sales & Marketing

Meril Healthcare Pvt. Ltd.
612, B-Wing, Bonanaza
Sahar Plaza, Andheri East
Mumbai 400 059
Maharashtra. India.
T: +91 22 3935 0700
F: +91 22 4047 9717

India Sales & Marketing

Meril Healthcare Pvt. Ltd.
512, Midas, Sahar Plaza, J.B. Nagar
Andheri East, Mumbai 400 059
Maharashtra. India.
T: +91 22 4047 9797

Meril South America

DOC MED LTDA
1079 – Cep: 04077-003 - Moema
Sao Paulo, Brazil
T/F: +55 11 3624 5935/6

Meril GmbH

Bornheimer Strasse 135 – 137
D – 53119 Bonn, Germany
T/F: +49 228 7100 4000/1

Meril Tibbi Cihazlar

Imalat Ve Ticaret A.S.
6, Mimar Sinan Mah.,
Cavusbasi Cad, Ozde Sok
Aydin Eksi Is Merkezi Kat:1
Cekmekoy/Istanbul, Turkey
T: +90 53 2272 5172



Health Care

CLAVO™ Catalogue

Nail | Screw | Instrument

Meril Orthopedic

Meril is known for an uncompromising attitude when it comes to quality of design. Providing the medical world with zero defect devices is our only motto. Our orthopedic trauma products---Meril Locking Plate System 3.5, Meril Locking Plate System 4.5, Meril Compression Plate System 3.5, Meril Compression Plate System 4.5, Meril Dynamic Hip Fixation System, and Meril Locking Distal Radial System 2.4 are specifically designed to treat fractures, bone tumour and deformities. They boast of high fatigue strength, proper fixation and negligible soft tissue damage.

We have been providing trauma implants for the foot, shoulder, hip, hand, pelvis, and condyles. In order to ensure smooth and trouble-free operating procedures leading medical professionals rely only on our Trauma products. We have carried out performance testing to prove our design for its endurance and durability and the outcome of these tests prove that they are going to create a dynamic image in the history of joint replacement.



Index

- CLAVO™ Proximal Femoral Rotational Stability Nail 1
- CLAVO™ Cannulated Titanium Tibial Nail 19
- CLAVO™ Elastic Titanium Nail 39
- CLAVO™ Shaft Femur Titanium Nail Cannulated 45
- CLAVO™ Shaft Tibia Titanium Nail Cannulated 46

CLAVO™ Proximal Femoral Rotational Stability Nail

Proximal Femoral Rotational - Stability Nail

Material : TAN(Ti-6Al-7Nb)

PRODUCT CODE	DIAMETER	LENGTH
--------------	----------	--------

MT-NT0109170	Ø 9.0	170
MT-NT0110170	Ø 10.0	170
MT-NT0111170	Ø 11.0	170
MT-NT0112170	Ø 12.0	170
MT-NT0109200	Ø 9.0	200
MT-NT0110200	Ø 10.0	200
MT-NT0111200	Ø 11.0	200
MT-NT0112200	Ø 12.0	200
MT-NT0109240	Ø 9.0	240
MT-NT0110240	Ø 10.0	240
MT-NT0111240	Ø 11.0	240
MT-NT0112240	Ø 12.0	240
MT-NT0109300L	Ø 9.0	300
MT-NT0109320L	Ø 9.0	320
MT-NT0109340L	Ø 9.0	340
MT-NT0109360L	Ø 9.0	360
MT-NT0109380L	Ø 9.0	380
MT-NT0109400L	Ø 9.0	400
MT-NT0109420L	Ø 9.0	420
MT-NT0109440L	Ø 9.0	440
MT-NT0110300L	Ø 10.0	300
MT-NT0110320L	Ø 10.0	320
MT-NT0110340L	Ø 10.0	340
MT-NT0110360L	Ø 10.0	360
MT-NT0110380L	Ø 10.0	380
MT-NT0110400L	Ø 10.0	400
MT-NT0110420L	Ø 10.0	420
MT-NT0110440L	Ø 10.0	440
MT-NT0112300L	Ø 12.0	300
MT-NT0112320L	Ø 12.0	320
MT-NT0112340L	Ø 12.0	340
MT-NT0112360L	Ø 12.0	360
MT-NT0112380L	Ø 12.0	380
MT-NT0112400L	Ø 12.0	400
MT-NT0112420L	Ø 12.0	420
MT-NT0112440L	Ø 12.0	440
MT-NT0109300R	Ø 9.0	300
MT-NT0109320R	Ø 9.0	320
MT-NT0109340R	Ø 9.0	340
MT-NT0109360R	Ø 9.0	360
MT-NT0109380R	Ø 9.0	380
MT-NT0109400R	Ø 9.0	400
MT-NT0109420R	Ø 9.0	420
MT-NT0109440R	Ø 9.0	440
MT-NT0110300R	Ø 10.0	300



CLAVO™ Proximal Femoral Rotational Stability Nail

Proximal Femoral Rotational - Stability Nail

Material : TAN(Ti-6Al-7Nb)

PRODUCT CODE	DIAMETER	LENGTH
--------------	----------	--------

MT-NT0110320R	Ø 10.0	320
MT-NT0110340R	Ø 10.0	340
MT-NT0110360R	Ø 10.0	360
MT-NT0110380R	Ø 10.0	380
MT-NT0110400R	Ø 10.0	400
MT-NT0110420R	Ø 10.0	420
MT-NT0110440R	Ø 10.0	440
MT-NT0112300R	Ø 12.0	300
MT-NT0112320R	Ø 12.0	320
MT-NT0112340R	Ø 12.0	340
MT-NT0112360R	Ø 12.0	360
MT-NT0112380R	Ø 12.0	380
MT-NT0112400R	Ø 12.0	400
MT-NT0112420R	Ø 12.0	420
MT-NT0112440R	Ø 12.0	440



CLAVO™ Proximal Femoral Rotational Stability Nail

Anti-Rotation Screw for PFRN

Material : TAN(Ti-6Al-7Nb)

PRODUCT CODE	DIAMETER	LENGTH
--------------	----------	--------

MT-BT0100075	12	75
MT-BT0100080	12	80
MT-BT0100085	12	85
MT-BT0100090	12	90
MT-BT0100095	12	95
MT-BT0100100	12	100
MT-BT0100105	12	105
MT-BT0100110	12	110
MT-BT0100115	12	115
MT-BT0100120	12	120



Nail Cap for PFRN

Material : TAV(Ti-6Al-4V)

PRODUCT CODE	DIAMETER	LENGTH
--------------	----------	--------

MT-ET0100000	12	0
MT-ET0100005	12	5
MT-ET0100010	12	10
MT-ET0100015	12	15




CLAVO™ Proximal Femoral Rotational Stability Nail

4.9mm Locking Bolt, Self Tapping

Material : TAN(Ti-6Al-7Nb) | Head Diameter : 8.0mm | Outer Diameter : 4.9mm
Core Diameter : 4.25mm | Pitch : 2.75mm

PRODUCT CODE	LENGTH
--------------	--------

MT-ST2549026	26
MT-ST2549028	28
MT-ST2549030	30
MT-ST2549032	32
MT-ST2549034	34
MT-ST2549036	36
MT-ST2549038	38
MT-ST2549040	40
MT-ST2549042	42
MT-ST2549044	44
MT-ST2549046	46
MT-ST2549048	48
MT-ST2549050	50
MT-ST2549052	52
MT-ST2549054	54
MT-ST2549056	56
MT-ST2549058	58
MT-ST2549060	60
MT-ST2549064	64
MT-ST2549068	68
MT-ST2549072	72
MT-ST2549076	76
MT-ST2549080	80
MT-ST2549085	85
MT-ST2549090	90
MT-ST2549095	95
MT-ST2549100	100



CLAVO™ Proximal Femoral Rotational Stability Nail Instrument Set



CLAVO™ Proximal Femoral Rotational Stability Nail Instrument Set

Part No:MT-INPFRN055

Description

Ø11.0mm Reamer for PFRN Anti-rotation Screw



Part No:MT-INPFRN032

Description

Gear Sets for PFRN



Part No:MT-INPFRN002

Description

Drill Bit Ø11.0mm for PFRN Anti-rotation Screw



Part No:MT-INPFRN047

Description

Direct Measuring Device for guide wire Ø3.2mm



Part No:MT-INPFRN005

Description

Trocar Ø3.2mm for PFRN Anti-rotation Screw



Part No:MT-INPFRN006

Description

Drill Sleeve Ø3.2mm/11.0 for PFRN Anti-rotation Screw



CLAVO™ Proximal Femoral Rotational Stability Nail Instrument Set

Part No:MT-INPFRN007

Description

Protection Sleeve 16.0mm/11.0 for PFRN Anti-rotation Screw



Part No:MT-INPFRN033

Description

Buttress/Compression Nut for PFRN Anti-rotation Screw



Part No:MT-INPFRN057

Description

Impactor for PFRN Anti-rotation Screw



Part No:MT-INPFRN059

Description

Combined Hammer(500gm) for PFRN



Part No:MT-INPFRN008

Description

Cannulated Connecting Screw for PFRN



CLAVO™ Proximal Femoral Rotational Stability Nail Instrument Set

Part No:MT-INPFRN050

Description

Insertion Handle for PFRN



Part No:MT-INPFRN009

Description

AWL for PFRN



Part No:MT-INPFRN053

Description

T-Generic Interface Handle for PFRN



Part No:MT-INPFRN010

Description

Ø3.2mm Guide Wire for PFRN Anti-rotation Screw



CLAVO™ Proximal Femoral Rotational Stability Nail Instrument Set

Part No:MT-INPFRN003

Description

Drill Bit Ø17.0mm Flexible Cannulated for PFRN



Part No:MT-INPFRN011

Description

Protection Sleeve 20.0mm/17.0 for PFRN



Part No:MT-INPFRN012

Description

Drill sleeve for PFRN



Part No:MT-INPFRN034

Description

Hexagonal Wrench with T Handle



CLAVO™ Proximal Femoral Rotational Stability Nail Instrument Set

Part No:MT-INPFRN035

Description

Connector for PFRN



Part No:MT-INPFRN036

Description

Ø4.5mm Pin Wrench 120mm Length



Part No:MT-INPFRN052

Description

Holding Sleeve, Large



Part No:MT-INPFRN013

Description

Aiming Jig For Anti-Rotation Wire



Part No:MT-INPFRN014

Description

Drill Sleeve Ø3.2mm/5.6 for PFRN



CLAVO™ Proximal Femoral Rotational Stability Nail Instrument Set

Part No:MT-INPFRN015

Description

Protection Sleeve 11.0mm/8.0 for PFRN



Part No:MT-INPFRN016

Description

Drill Sleeve Ø4.0mm/8.0 for PFRN, Green



Part No:MT-INPFRN017

Description

Trocar Ø4.0mm for PFRN



Part No:MT-INPFRN048

Description

Measuring Device for Locking Bolt 20-100



Part No:MT-INPFRN049

Description

Radiographic Ruler for PFRN



Part No:MT-INPFRN018

Description

Ø2.8mm Guide Wire with Hook 460mm Length



CLAVO™ Proximal Femoral Rotational Stability Nail Instrument Set

Part No:MT-INPFRN001

Description

Drill Bit Ø4.0mm for PFRN Anti-rotation Screw



Part No:MT-INPFRN037

Description

Hexagonal Screw Driver Cannulated Ø6.5 & 7.3mm



Part No:MT-INPFRN038

Description

Hexagonal Screw Driver Large Ø3.5 & 300mm Length



Part No:MT-INPFRN019

Description

Static Lock Aiming Arm for Small & Extra Small PFRN Nail



Part No:MT-INPFRN054

Description

Screw Holding Forceps with self Holding for PFRN



CLAVO™ Proximal Femoral Rotational Stability Nail Instrument Set

Part No:MT-INPFRN020

Description

Dynamic Locking Aiming Arm for PFRN



Part No:MT-INPFRN039

Description

Hexagonal Spherical Wrench



Part No:MT-INPFRN040

Description

Soft Diamond for PFRN



Part No:MT-INPFRN041

Description

M8 Guide Rod for PFRN



Part No:MT-INPFRN042

Description

Limit Wrench for PFRN



Part No:MT-INPFRN022

Description

Guide Pin Ø2.5mmx450mm



CLAVO™ Proximal Femoral Rotational Stability Nail Instrument Set

Part No:MT-INPFRN043

Description

Clip - Fastner



Part No:MT-INPFRN023

Description

Sleeve for Positioning Rod



Part No:MT-INPFRN024

Description

Aiming Arm for Connecting Blocks - 260 to 340



Part No:MT-INPFRN025

Description

Aiming Arm for Connecting Blocks - 360 to 420



CLAVO™ Proximal Femoral Rotational Stability Nail Instrument Set

Part No:MT-INPFRN026

Description

Guide Rod for PFRN



Part No:MT-INPFRN027

Description

Remote Aiming Rack



Part No:MT-INPFRN028

Description

M5 Guide Rod Lock Wheel



Part No:MT-INPFRN029

Description

M6 Guide Rod Lock Wheel



CLAVO™ Proximal Femoral Rotational Stability Nail Instrument Set

Part No:MT-INPFRN044

Description

Ø9.0mm Positioning Rod



Part No:MT-INPFRN045

Description

Ø10.0mm Positioning Rod



Part No:MT-INPFRN030

Description

Drill Sleeve 8.0/6.0



Part No:MT-INPFRN031

Description

Ø6.0mm Trocar



Part No:MT-INPFRN056

Description

Ø6.0mm Bone Reamer



CLAVO™ Proximal Femoral Rotational Stability Nail Instrument Set

Part No:MT-INPFRN004

Description

Ø6.0mm Bone Drill



Part No:MT-INPFRN046

Description

Allen Wrench



Part No:MT-INPFRN060

Description


Wrench for PFRN



Part No:MT-INPFRN058

Description


Apparatus Box - Case



CLAVO™ Cannulated Titanium Tibial Nail

Cannulated Titanium Tibial Nail System | Material : TAN(Ti-6Al-7Nb)


PRODUCT CODE	DIAMETER	LENGTH
MT-NT0208255	Ø 8.0	255
MT-NT0208270	Ø 8.0	270
MT-NT0208285	Ø 8.0	285
MT-NT0208300	Ø 8.0	300
MT-NT0208315	Ø 8.0	315
MT-NT0208330	Ø 8.0	330
MT-NT0208345	Ø 8.0	345
MT-NT0208360	Ø 8.0	360
MT-NT0208375	Ø 8.0	375
MT-NT0208390	Ø 8.0	390
MT-NT0208405	Ø 8.0	405
MT-NT0208420	Ø 8.0	420
MT-NT0208435	Ø 8.0	435
MT-NT0208450	Ø 8.0	450
MT-NT0209255	Ø 9.0	255
MT-NT0209270	Ø 9.0	270
MT-NT0209285	Ø 9.0	285
MT-NT0209300	Ø 9.0	300
MT-NT0209315	Ø 9.0	315
MT-NT0209330	Ø 9.0	330
MT-NT0209345	Ø 9.0	345
MT-NT0209360	Ø 9.0	360
MT-NT0209375	Ø 9.0	375
MT-NT0209390	Ø 9.0	390
MT-NT0209405	Ø 9.0	405
MT-NT0209420	Ø 9.0	420
MT-NT0209435	Ø 9.0	435
MT-NT0209450	Ø 9.0	450
MT-NT0210255	Ø 10.0	255
MT-NT0210270	Ø 10.0	270
MT-NT0210285	Ø 10.0	285
MT-NT0210300	Ø 10.0	300
MT-NT0210315	Ø 10.0	315
MT-NT0210330	Ø 10.0	330
MT-NT0210345	Ø 10.0	345
MT-NT0210360	Ø 10.0	360
MT-NT0210375	Ø 10.0	375
MT-NT0210390	Ø 10.0	390
MT-NT0210405	Ø 10.0	405
MT-NT0210420	Ø 10.0	420
MT-NT0210435	Ø 10.0	435
MT-NT0210450	Ø 10.0	450
MT-NT0211255	Ø 11.0	255



CLAVO™ Cannulated Titanium Tibial Nail

Cannulated Titanium Tibial Nail System | Material : TAN(Ti-6Al-7Nb)

PRODUCT CODE	DIAMETER	LENGTH
MT-NT0211270	Ø 11.0	270
MT-NT0211285	Ø 11.0	285
MT-NT0211300	Ø 11.0	300
MT-NT0211315	Ø 11.0	315
MT-NT0211330	Ø 11.0	330
MT-NT0211345	Ø 11.0	345
MT-NT0211360	Ø 11.0	360
MT-NT0211375	Ø 11.0	375
MT-NT0211390	Ø 11.0	390
MT-NT0211405	Ø 11.0	405
MT-NT0211420	Ø 11.0	420
MT-NT0211435	Ø 11.0	435
MT-NT0211450	Ø 11.0	450
MT-NT0212255	Ø 12.0	255
MT-NT0212270	Ø 12.0	270
MT-NT0212285	Ø 12.0	285
MT-NT0212300	Ø 12.0	300
MT-NT0212315	Ø 12.0	315
MT-NT0212330	Ø 12.0	330
MT-NT0212345	Ø 12.0	345
MT-NT0212360	Ø 12.0	360
MT-NT0212375	Ø 12.0	375
MT-NT0212390	Ø 12.0	390
MT-NT0212405	Ø 12.0	405
MT-NT0212420	Ø 12.0	420
MT-NT0212435	Ø 12.0	435
MT-NT0212450	Ø 12.0	450
MT-NT0213255	Ø 13.0	255
MT-NT0213270	Ø 13.0	270
MT-NT0213285	Ø 13.0	285
MT-NT0213300	Ø 13.0	300
MT-NT0213315	Ø 13.0	315
MT-NT0213330	Ø 13.0	330
MT-NT0213345	Ø 13.0	345
MT-NT0213360	Ø 13.0	360
MT-NT0213375	Ø 13.0	375
MT-NT0213390	Ø 13.0	390
MT-NT0213405	Ø 13.0	405
MT-NT0213420	Ø 13.0	420
MT-NT0213435	Ø 13.0	435
MT-NT0213450	Ø 13.0	450



CLAVO™ Cannulated Titanium Tibial Nail

Nail Cap For CTTN | Material : TAN(Ti-6Al-7Nb)

PRODUCT CODE	DIAMETER	LENGTH
MT-ET0208000	-	-
MT-ET0208005	-	-
MT-ET0208010	-	-
MT-ET0208015	-	-
MT-ET0208000L	-	-
MT-ET0208005L	-	-
MT-ET0208010L	-	-
MT-ET0208015L	-	-



CLAVO™ Cannulated Titanium Tibial Nail

4.0mm Locking Screw Medullary Nail

Material : TAN(Ti-6Al-7Nb) | Head Diameter : 8.0mm | Outer Diameter : 3.9mm
Core Diameter : 3.25mm | Pitch : 1.25mm

PRODUCT CODE	LENGTH
MT-ST1540018	18mm
MT-ST1540020	20mm
MT-ST1540022	22mm
MT-ST1540024	24mm
MT-ST1540026	26mm
MT-ST1540028	28mm
MT-ST1540030	30mm
MT-ST1540032	32mm
MT-ST1540034	34mm
MT-ST1540036	36mm
MT-ST1540038	38mm
MT-ST1540040	40mm
MT-ST1540042	42mm
MT-ST1540044	44mm
MT-ST1540046	46mm
MT-ST1540048	48mm
MT-ST1540050	50mm
MT-ST1540052	52mm
MT-ST1540054	54mm
MT-ST1540056	56mm
MT-ST1540058	58mm
MT-ST1540060	60mm
MT-ST1540062	62mm
MT-ST1540064	64mm
MT-ST1540066	66mm
MT-ST1540068	68mm
MT-ST1540070	70mm
MT-ST1540072	72mm
MT-ST1540074	74mm
MT-ST1540076	76mm
MT-ST1540078	78mm
MT-ST1540080	80mm



CLAVO™ Cannulated Titanium Tibial Nail

5.0mm Locking Screw Medullary Nail


Material : TAN(Ti-6Al-7Nb) | Head Diameter : 8.0mm | Outer Diameter : 4.9mm
Core Diameter : 4.25mm | Pitch : 1.25mm

PRODUCT CODE	LENGTH
--------------	--------

MT-ST2750026	26mm
MT-ST2750028	28mm
MT-ST2750030	30mm
MT-ST2750032	32mm
MT-ST2750034	34mm
MT-ST2750036	36mm
MT-ST2750038	38mm
MT-ST2750040	40mm
MT-ST2750042	42mm
MT-ST2750044	44mm
MT-ST2750046	46mm
MT-ST2750048	48mm
MT-ST2750050	50mm
MT-ST2750052	52mm
MT-ST2750054	54mm
MT-ST2750056	56mm
MT-ST2750058	58mm
MT-ST2750060	60mm
MT-ST2750062	62mm
MT-ST2750064	64mm
MT-ST2750066	66mm
MT-ST2750068	68mm
MT-ST2750070	70mm
MT-ST2750072	72mm
MT-ST2750074	74mm
MT-ST2750076	76mm
MT-ST2750078	78mm
MT-ST2750080	80mm
MT-ST2750085	85mm
MT-ST2750090	90mm
MT-ST2750095	95mm
MT-ST2750100	100mm



CLAVO™ Cannulated Titanium Tibial Nail Instrument Set



CLAVO™ Cannulated Titanium Tibial Nail Instrument Set

Part No:MT-INETNS001

Description

Radiographic Ruler



Part No:MT-INETNS002

Description

Guide wire with threaded tip, Ø3.2, length 320mm



Part No:MT-INETNS003

Description

Holder for guide Wire Reaming Drill, Proximal



Part No:MT-INETNS004

Description

Protection Sleeve



Part No:MT-INETNS005

Description

Chisel Cannulated



CLAVO™ Cannulated Titanium Tibial Nail Instrument Set

Part No:MT-INETNS006

Description

Tissue Protector



Part No:MT-INETNS007

Description

Awl, Cannulated



Part No:MT-INETNS008

Description

Reaming Drill, Proximal



Part No:MT-INETNS009

Description

T-Handle with quick coupling



CLAVO™ Cannulated Titanium Tibial Nail Instrument Set

Part No:MT-INETNS010

Description

Reduction Rod



Part No:MT-INETNS011

Description

Guide wire with olive head, Ø2.5 length 1000mm



Part No:MT-INETNS012

Description

Insertion Handle



Part No:MT-INETNS013

Description

Connecting Screw for insertion handle



CLAVO™ Cannulated Titanium Tibial Nail Instrument Set

Part No:MT-INETNS014

Description

Wrench for connecting Screw



Part No:MT-INETNS015

Description

Connector Block



Part No:MT-INETNS016

Description

Connector Shaft



Part No:MT-INETNS017

Description

Ratchet Wrench, Sw11



CLAVO™ Cannulated Titanium Tibial Nail Instrument Set

Part No:MT-INETNS018

Description

Extraction Screw for Tibial Nail



Part No:MT-INETNS019

Description

Slide Hammer



Part No:MT-INETNS020

Description

Aiming Arm



Part No:MT-INETNS021

Description


Connecting Screw, for aiming shaft



Part No:MT-INETNS022

Description

Wrench for connecting Screw, Sw5



CLAVO™ Cannulated Titanium Tibial Nail Instrument Set

Part No:MT-INETNS023

Description

Drill Sleeve for Calibrating Pin, Ø5.0



Part No:MT-INETNS024

Description


Drill Bit Ø5.0, length 250mm



Part No:MT-INETNS025

Description


Drill Bit with flat head



Part No:MT-INETNS026

Description

Calibrating Pin



CLAVO™ Cannulated Titanium Tibial Nail Instrument Set

Part No:MT-INETNS027

Description

Stretch Aiming Arm, Distal



Part No:MT-INETNS028

Description

Aiming Arm Distal



Part No:MT-INETNS029

Description

Connecting Screw for Aiming Shaft, M6



Part No:MT-INETNS030

Description

Protection Sleeve, Ø8.1/ Ø10




CLAVO™ Cannulated Titanium Tibial Nail Instrument Set

Part No:MT-INETNS031

Description


Trocar Ø8.1



Part No:MT-INETNS032

Description

Dril Sleeve for Drill Bits, Ø8.1/ Ø4.2



Part No:MT-INETNS033

Description

Drill Bit, Ø4.2 Length 300mm



Part No:MT-INETNS034

Description

Drill Bit, Ø4.2 Length 350mm



Part No:MT-INETNS035

Description

Fixation Sleeve, Ø4.2




CLAVO™ Cannulated Titanium Tibial Nail Instrument Set

Part No:MT-INETNS036

Description

Wrench for Fixation Sleeve, Sw3



Part No:MT-INETNS037

Description

Depth Gauge



Part No:MT-INETNS038

Description

Screwdriver, Stardrive, T25



Part No:MT-INETNS039

Description

Screwdriver, Stardrive, T25



Part No:MT-INETNS040

Description

Proximal Aiming Arm, Left



CLAVO™ Cannulated Titanium Tibial Nail Instrument Set

Part No:MT-INETNS041

Description

Proximal Aiming Arm, Right



Part No:MT-INETNS042

Description

Aiming Arm for Cancellous Bone, Proximal



Part No:MT-INETNS043

Description

Drill Sleeve, Ø8.1/ Ø3.4



Part No:MT-INETNS044

Description

Fixation Sleeve, Ø3.4




CLAVO™ Cannulated Titanium Tibial Nail Instrument Set

Part No:MT-INETNS045

Description

Compression Screw



Part No:MT-INETNS046

Description

Screwdriver , stardrive, cannulated T25



Part No:MT-INETNS047

Description


Universal Wrench, SW6.5



Part No:MT-INETNS048

Description

Screwdriver Shaft for End Cap



Part No:MT-INETNS049

Description

U - Clip



CLAVO™ Cannulated Titanium Tibial Nail Instrument Set

Part No:MT-INETNS050

Description

Drill Bit Ø3.4, length 280mm



Part No:MT-INETNS051

Description


Drill Bit, Ø3.4 length 320mm



Part No:MT-INETNS052

Description

Instrument Case



CLAVO™ Elastic Titanium Nail

Elastic Titanium Nail | Material : TAV(Ti-6Al-4V)

PRODUCT CODE	DIAMETER	LENGTH
MT-NT0515300	Ø1.5	300
MT-NT0520440	Ø2.0	440
MT-NT0525440	Ø2.5	440
MT-NT0530440	Ø3.0	440
MT-NT0535440	Ø3.5	440
MT-NT0540440	Ø4.0	440



CLAVO™ Elastic Titanium Nail Instrument Set



CLAVO™ Elastic Titanium Nail Instrument Set

Part No:MT-INTELN001

Description

Awl for TEN



Part No:MT-INTELN002

Description

Insertion for TEN



Part No:MT-INTELN003

Description

Impactor for TEN



Part No:MT-INTELN004

Description

Hammer for TEN




CLAVO™ Elastic Titanium Nail Instrument Set

Part No:MT-INTELN005

Description

Extraction Pliers for TEN



Part No:MT-INTELN006

Description

Drill Bit Ø2.5 length 127mm



Part No:MT-INTELN007

Description

Drill Bit Ø3.2 length 150mm



Part No:MT-INTELN008

Description

Drill Bit Ø4.5, length 147mm



Part No:MT-INLFC001

Description

Double Drill Guide, Ø3.2/ Ø4.5



CLAVO™ Elastic Titanium Nail Instrument Set

Part No:MT-INTELN010

Description

Screwdriver Shaft for end cap, Ø2.0/ Ø2.5



Part No:MT-INTELN011

Description

Screwdriver Shaft for end cap, Ø3.0/ Ø3.5/ Ø4.5



Part No:MT-INTELN013

Description

Pin Wrench



Part No:MT-INTELN012

Description

Instrument Case



CLAVO™ Shaft Femur Titanium Nail Cannulated

Shaft Femur Titanium Nail Cannulated | Material : TAN(ASTM-F1295-11)


PRODUCT CODE	DIAMETER	LENGTH
MT-NT1309340	Ø 9.0	340mm
MT-NT1309360	Ø 9.0	360mm
MT-NT1309380	Ø 9.0	380mm
MT-NT1309400	Ø 9.0	400mm
MT-NT1309420	Ø 9.0	420mm
MT-NT1309440	Ø 9.0	440mm
MT-NT1310340	Ø 10.0	340mm
MT-NT1310360	Ø 10.0	360mm
MT-NT1310380	Ø 10.0	380mm
MT-NT1310400	Ø 10.0	400mm
MT-NT1310420	Ø 10.0	420mm
MT-NT1310440	Ø 10.0	440mm
MT-NT1311340	Ø 11.0	340mm
MT-NT1311360	Ø 11.0	360mm
MT-NT1311380	Ø 11.0	380mm
MT-NT1311400	Ø 11.0	400mm
MT-NT1311420	Ø 11.0	420mm
MT-NT1311440	Ø 11.0	440mm
MT-NT1312340	Ø 12.0	340mm
MT-NT1312360	Ø 12.0	360mm
MT-NT1312380	Ø 12.0	380mm
MT-NT1312400	Ø 12.0	400mm
MT-NT1312420	Ø 12.0	420mm
MT-NT1312440	Ø 12.0	440mm



CLAVO™ Shaft Tibia Titanium Nail Cannulated

Shaft Tibia Titanium Nail Cannulated | Material : TAN(ASTM-F1295-11)

PRODUCT CODE	DIAMETER	LENGTH
MT-NT1408280	Ø 8.0	280mm
MT-NT1408300	Ø 8.0	300mm
MT-NT1408320	Ø 8.0	320mm
MT-NT1408340	Ø 8.0	340mm
MT-NT1408360	Ø 8.0	360mm
MT-NT1408380	Ø 8.0	380mm
MT-NT1409280	Ø 9.0	280mm
MT-NT1409300	Ø 9.0	300mm
MT-NT1409320	Ø 9.0	320mm
MT-NT1409340	Ø 9.0	340mm
MT-NT1409360	Ø 9.0	360mm
MT-NT1409380	Ø 9.0	380mm
MT-NT1410280	Ø 10.0	280mm
MT-NT1410300	Ø 10.0	300mm
MT-NT1410320	Ø 10.0	320mm
MT-NT1410340	Ø 10.0	340mm
MT-NT1410360	Ø 10.0	360mm
MT-NT1410380	Ø 10.0	380mm



3.9mm Locking Bolt Self Tapping

Material : TAV(Ti-6Al-4V) | Head Diameter : 8.0mm | Outer Diameter : 3.9mm
 Core Diameter : 3.25mm | Pitch : 2.3mm

PRODUCT CODE	LENGTH
MT-ST1339022	22
MT-ST1339024	24
MT-ST1339026	26
MT-ST1339028	28
MT-ST1339030	30
MT-ST1339032	32
MT-ST1339034	34
MT-ST1339036	36
MT-ST1339038	38
MT-ST1339040	40
MT-ST1339042	42
MT-ST1339044	44
MT-ST1339046	46
MT-ST1339048	48
MT-ST1339050	50
MT-ST1339052	52
MT-ST1339054	54
MT-ST1339056	56
MT-ST1339058	58
MT-ST1339060	60



4.9mm Locking Bolt, Self Tapping

Material : TAN(Ti-6Al-7Nb) | Head Diameter : 8.0mm | Outer Diameter : 4.9mm
 Core Diameter : 4.25mm | Pitch : 2.75mm

PRODUCT CODE	LENGTH
MT-ST2549026	26
MT-ST2549028	28
MT-ST2549030	30
MT-ST2549032	32
MT-ST2549034	34
MT-ST2549036	36
MT-ST2549038	38
MT-ST2549040	40
MT-ST2549042	42
MT-ST2549044	44
MT-ST2549046	46
MT-ST2549048	48
MT-ST2549050	50
MT-ST2549052	52
MT-ST2549054	54
MT-ST2549056	56
MT-ST2549058	58
MT-ST2549060	60
MT-ST2549064	64
MT-ST2549068	68
MT-ST2549072	72
MT-ST2549076	76
MT-ST2549080	80
MT-ST2549085	85
MT-ST2549090	90
MT-ST2549095	95
MT-ST2549100	100mm

