


MIRUSTM

Linear Cutter



MIRUS™

Linear Cutter

Cut Line Marking

Cartridge Half

Anvil Tip

Anvil Half

Intermediate Alignment



- Uniform compression & consistent staple formation
- Low profile anvil tip to minimise tissue trauma



Cut-line & Staple Line markers for precise cutting & stapling




New knife with every reload guarantees smooth & precise transection

Firm Ergonomic Grip

Dual Side Firing Knob

Push Button

Rear Hinge



Every reload has a safety lock to prevent firing over empty reload



- Push-Button for quick release
- Ambidextrous-Dual side firing for ease of use



Rear hinge alignment offers ease of handling

Technical Specifications

Stapler	Description	Relods	Open Length	Close Length	Color	Number of Firings	Number of Staples	Rows of Staple	Staple Line Length	Cut Length
MLC60	Mirus Disp. Linear Cutter 60	MLCR60-3.8 MLCR60-4.8	3.8 mm 4.8 mm	1.5 mm 2.0 mm	Blue Green	8	64	4	64 mm	60 mm
MLC80	Mirus Disp. Linear Cutter 80	MLCR80-3.8 MLCR80-4.8	3.8 mm 4.8 mm	1.5 mm 2.0 mm	Blue Green	8	84	4	84 mm	80 mm
MLC100	Mirus Disp. Linear Cutter 100	MLCR100-3.8 MLCR100-4.8	3.8 mm 4.8 mm	1.5 mm 2.0 mm	Blue Green	8	104	4	104 mm	100 mm

Meril's Global Presence

Meril Tibbi Cihazlar

İmalat ve Ticaret A.Ş.
İçerenköy Mah. Çetinkaya Sok.
Prestij Plaza No:28, Kat:4 Ataşehir,
34752 İstanbul / Turkey
T : +90 216 641 44 24
F : +90 216 641 44 25

Meril GmbH.

Bornheimer Strasse 135-137,
D-53119 Bonn.
Germany.
T : +49 228 7100 4000
F : +49 228 7100 4001

Meril, inc.

2436 Emrick Boulevard,
Bethlehem, PA - 18020
T : (610)500-2080
F : (610)317-1672

Meril South America

Doc Med LTDA Al. dos Tupiniquins,
1079 - Cep: 04077-003 - Moema.
Sao Paulo. Brazil.
T : +55 11 3624 5935
F : +55 11 3624 5936

Meril Medical LLC.

Nauchnyi Proezd 19,
Moscow , Russia – 117 246.
Office - +7 495 772 7643

Meril UK Subsidiary

Crowe UK LLP, St. Brides
House, 10 Salisbury Square,
London EC4Y 8EH, UK

EU Representative.

Obelis S.A.Bd, General Wahis 53,
1030, Brussels, Belgium.
T : +32 2 732 5954
F : +32 2 732 6003
E : mail@obelis.net

Meril SA Pty Ltd

102, 104, S101 and S102, Boulevard
West Office Park, 142 Western
Service Road, Erf 813
Woodmead Extension 17 Sandton,
Johannesburg – 2191 South Africa
T : +27 11 465-2049
F : +27 86 471 7941

Meril Endo Surgery Pvt. Ltd

Third Floor E1-E3, Meril Park,
Survey Number 135/2/B & 172/2
Muktanand Marg, Chala,
Vapi 396 191, Gujarat. India.
T : +91 260 305 2100
F : +91 260 305 2125
E : askinfo@merillife.com

Meril China Co. Ltd.

2301b,23f, Lixin Plaza,
no 90, South Hubin Road,
Xiamen, China.
T : 0086-592-5368505
F : 0086-592-5368519

Meril Bangladesh Pvt. Ltd.

Union Heights-01,
55-2, Bir Uttam Qazi
Nuruzzaman Sarak,
West Panthapath,
Dhaka-1205, Bangladesh
T : + 88 02 9122714/24

Meril Life Sciences Pvt. Ltd.

301, A - Wing, Business Square,
Chakala, Andheri Kurla Road,
Andheri East, Mumbai 400 093
T : + 91 22 39350700
F : + 91 22 39350777
E : askinfo@merillife.com

E askinfo@merillife.com
W www.merillife.com

Meril

Meril Endo-Surgery Pvt. Ltd.
E1-E3, Meril Park, Survey No. 135/2/B & 174/2, Muktanand Marg, Chala,
Vapi - 396 191, Gujarat, India